

β -Tubulin Rabbit pAb

Catalog No.: A1010

Basic Information

Catalog No.

A1010

Observed MW

55kDa

Calculated MW

50kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Rat

Recommended Dilutions

WB 1:500 - 1:2000

Background

This gene encodes a beta tubulin protein. This protein forms a dimer with alpha tubulin and acts as a structural component of microtubules. Mutations in this gene cause cortical dysplasia, complex, with other brain malformations 6. Alternative splicing results in multiple splice variants. There are multiple pseudogenes for this gene on chromosomes 1, 6, 7, 8, 9, and 13.

Immunogen Information

Gene ID

203068

Swiss Prot

P07437


Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 40-260 of human β -Tubulin (NP_006077.2).

Synonyms

CDCBM6; CSCSC1; M40; OK/SW-cl.56; TUBB1; TUBB5; Beta Tubulin; TUBB; β -Tubulin

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.