Leader in Biomolecular Solutions for Life Science

Beclin 1 Mouse mAb

Catalog No.: A10101 2 Publications



Basic Information

Observed MW

60kDa

Calculated MW

52kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human, Mouse, Rat, Hamster, Monkey

Background

This gene encodes a protein that regulates autophagy, a catabolic process of degradation induced by starvation. The encoded protein is a component of the phosphatidylinositol-3-kinase (PI3K) complex which mediates vesicle-trafficking processes. This protein is thought to play a role in multiple cellular processes, including tumorigenesis, neurodegeneration and apoptosis. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB

1:200 - 1:1000

Immunogen Information

Gene ID

Swiss Prot Q14457

8678

Immunogen

Recombinant protein of human Beclin 1

Synonyms

ATG6; VPS30; beclin1; Beclin 1

Contact

•

www.abclonal.com

Product Information


Source Mouse **Isotype** IgG **Purification**Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using Beclin 1 Mouse mAb (A10101) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS003) at 1:10000 dilution.

Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021).

Exposure time: 90s.