

RPRD1B Rabbit pAb

Catalog No.: A10102

Basic Information

Catalog No.

A10102

Observed MW

37kDa

Calculated MW

36kDa

Category

Primary antibody

Applications

WB, IHC

Cross-Reactivity

Human, Mouse, Rat

Recommended Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200

Background

Immunogen Information

Gene ID

58490

Swiss Prot

Q9NQG5

Immunogen

Recombinant protein of human rprd1b

Synonyms

RPRD1B;C20orf77;CREPT;NET60;dj1057B20.2;K-H;Kub5-Hera;rprd1b

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using rprd1b antibody (A10102) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 15s.