

GZMK Rabbit pAb

Catalog No.: A10219

Basic Information

Catalog No.

A10219

Observed MW

29kDa

Calculated MW

28kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Mouse, Rat

Recommended Dilutions

WB 1:500 - 1:2000

Background

This gene product is a member of a group of related serine proteases from the cytoplasmic granules of cytotoxic lymphocytes. Cytolytic T lymphocytes (CTL) and natural killer (NK) cells share the remarkable ability to recognize, bind, and lyse specific target cells. They are thought to protect their host by lysing cells bearing on their surface 'nonself' antigens, usually peptides or proteins resulting from infection by intracellular pathogens. The protein described here lacks consensus sequences for N-glycosylation present in other granzymes.

Immunogen Information

Gene ID

3003

Swiss Prot

P49863

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 27-264 of human GZMK (NP_002095.1).

Synonyms

GZMK;TRYP2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using GZMK antibody (A10219) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 10s.