

A10678

Leader in Biomolecular Solutions for Life Science



AHI1 Rabbit pAb

Catalog No.: A10678

Basic Information

Observed MW

Calculated MW
137kDa

Category
Polyclonal Antibody

Applications
ELISA

Cross-Reactivity
Mouse,Rat

Background

This gene is apparently required for both cerebellar and cortical development in humans. This gene mutations cause specific forms of Joubert syndrome-related disorders. Joubert syndrome (JS) is a recessively inherited developmental brain disorder with several identified causative chromosomal loci. Alternatively spliced transcript variants encoding different isoforms have been identified.

Recommended Dilutions

Immunogen Information

Gene ID	Swiss Prot
54806	Q8N157

Immunogen
Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human AHI1 (NP_001128304.1).

Synonyms
ORF1; AHI-1; JBTS3; dj71N10.1; AHI1

Contact

 | www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.