

A11957

Leader in Biomolecular Solutions for Life Science



ADGRB2 Rabbit pAb

Catalog No.: A11957

Basic Information

Observed MW

Refer to Figures

Calculated MW

173kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene encodes a seven-span transmembrane protein that is thought to be a member of the secretin receptor family. The encoded protein is a brain-specific inhibitor of angiogenesis. The mature peptide may be further cleaved into additional products (PMID:20367554). Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

576

Swiss Prot

O60241

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 700-800 of human ADGRB2 (NP_001281264.1).

Synonyms

BAI2; ADGRB2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.