

A12127

Leader in Biomolecular Solutions for Life Science



# ACAP2 Rabbit pAb

Catalog No.: A12127

## Basic Information

---

### Observed MW

115kDa

### Calculated MW

88kDa

### Category

Polyclonal Antibody

### Applications

WB,IF/ICC,ELISA

### Cross-Reactivity

Human,Mouse

## Background

---

Enables GTPase activator activity. Acts upstream of or within actin filament-based process. Located in ruffle.

## Recommended Dilutions

---

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:200

## Immunogen Information

---

### Gene ID

23527

### Swiss Prot

Q15057

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 50-160 of human ACAP2 (NP\_036419.3).

### Synonyms

CENTB2; CNT-B2; ACAP2

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

## Validation Data



Western blot analysis of various lysates using ACAP2 Rabbit pAb (A12127) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 5s.



Immunofluorescence analysis of U2OS cells using ACAP2 Rabbit pAb (A12127) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.