

A13183

Leader in Biomolecular Solutions for Life Science



## SAAL1 Rabbit pAb

Catalog No.: A13183

### Basic Information

**Observed MW**

54-58kDa

**Calculated MW**

54kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human, Mouse, Rat

### Background

Acts upstream of or within positive regulation of synoviocyte proliferation. Located in nucleoplasm.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

113174

**Swiss Prot**

Q96ER3

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 175-474 of human SAAL1 (NP\_612430.2).

**Synonyms**

SPACIA1; SAAL1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**


Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using SAAL1 Rabbit pAb (A13183) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Enhanced Kit (RM00021).  
Exposure time: 10s.