

# TESC Rabbit pAb

Catalog No.: A13219

## Basic Information

### Observed MW

Refer to Figures

### Calculated MW

25kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse

## Background

Enables calcium ion binding activity. Involved in several processes, including cellular response to retinoic acid; positive regulation of macromolecule metabolic process; and positive regulation of myeloid cell differentiation. Located in several cellular components, including cytosol; lamellipodium; and nucleoplasm.

## Recommended Dilutions

WB 1:1000 - 1:2000

## Immunogen Information

### Gene ID

54997

### Swiss Prot

Q96BS2

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 135-214 of human TESC (NP\_060369.3).

### Synonyms

TSC; CHP3; TESC

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.