

A13255

Leader in Biomolecular Solutions for Life Science



## USP6 Rabbit pAb

Catalog No.: A13255

### Basic Information

#### Observed MW

Refer to Figures

#### Calculated MW

159kDa

#### Category

Mouse Monoclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

Enables thiol-dependent deubiquitinase. Involved in protein deubiquitination and regulation of vesicle-mediated transport. Located in plasma membrane and recycling endosome.

### Recommended Dilutions

WB 1:1000 - 1:2000

### Immunogen Information

#### Gene ID

9098

#### Swiss Prot

P35125

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 55-120 of human USP6 (NP\_004496.2).

#### Synonyms

HRP1; TRE2; TRE17; Tre-2; TRESMCR; USP6-short; USP6

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.