

USP6 Rabbit pAb

Catalog No.: A13255

Basic Information

Observed MW

Refer to Figures

Calculated MW

159kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Enables thiol-dependent deubiquitinase. Involved in protein deubiquitination and regulation of vesicle-mediated transport. Located in plasma membrane and recycling endosome.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

9098

Swiss Prot

P35125

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 55-120 of human USP6 (NP_004496.2).

Synonyms

HRP1; TRE2; TRE17; Tre-2; TRESMCR; USP6-short; USP6

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.