

A13264

Leader in Biomolecular Solutions for Life Science



# AGER Rabbit pAb

Catalog No.: A13264 **7 Publications**

## Basic Information

### Observed MW

45kDa/55kDa

### Calculated MW

43kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Mouse,Rat

## Background

The advanced glycosylation end product (AGE) receptor encoded by this gene is a member of the immunoglobulin superfamily of cell surface receptors. It is a multiligand receptor, and besides AGE, interacts with other molecules implicated in homeostasis, development, and inflammation, and certain diseases, such as diabetes and Alzheimer's disease. Many alternatively spliced transcript variants encoding different isoforms, as well as non-protein-coding variants, have been described for this gene (PMID:18089847).

## Recommended Dilutions

WB 1:5000 - 1:10000

## Immunogen Information

### Gene ID

177

### Swiss Prot

Q15109

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 40-210 of human AGER (NP\_001127.1).

### Synonyms

RAGE; sRAGE; SCARJ1; AGER

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of lysates from Rat lung, using AGER Rabbit pAb (A13264) at 1:5000/1:10000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.