

A13408

Leader in Biomolecular Solutions for Life Science



SMPX Rabbit pAb

Catalog No.: A13408

Basic Information

Observed MW

Calculated MW
10kDa

Category
Polyclonal Antibody

Applications
IF/ICC,ELISA

Cross-Reactivity
Human

Background

This gene encodes a small protein that has no known functional domains. Mutations in this gene are a cause of X-linked deafness-4, and the encoded protein may play a role in the maintenance of inner ear cells subjected to mechanical stress. Alternatively spliced transcript variants have been observed for this gene.

Recommended Dilutions

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID	Swiss Prot
23676	Q9UHP9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-88 of human SMPX (NP_055147.1).

Synonyms

Csl; DFN6; MPD7; DFNX4; Chisel; SMPX

Contact



www.abclonal.com


Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of HeLa cells using SMPX antibody (A13408). Blue: DAPI for nuclear staining.