

A13637

Leader in Biomolecular Solutions for Life Science



TNFAIP3 Rabbit pAb

Catalog No.: A13637 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

90kDa

Category

Polyclonal Antibody

Applications

WB,IF/ICC,ELISA

Cross-Reactivity

Human,Rat

Background

This gene was identified as a gene whose expression is rapidly induced by the tumor necrosis factor (TNF). The protein encoded by this gene is a zinc finger protein and ubiquitin-editing enzyme, and has been shown to inhibit NF-kappa B activation as well as TNF-mediated apoptosis. The encoded protein, which has both ubiquitin ligase and deubiquitinase activities, is involved in the cytokine-mediated immune and inflammatory responses. Several transcript variants encoding the same protein have been found for this gene.

Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200

Immunogen Information

Gene ID	Swiss Prot
7128	P21580

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human TNFAIP3 (NP_006281.1).

Synonyms

A20; AISBL; AIFBL1; OTUD7C; TNFA1P2; TNFAIP3

Contact



www.abclonal.com

Product Information


Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of U2OS cells using TNFAIP3 antibody (A13637). Blue: DAPI for nuclear staining.