

KIAA1456 Rabbit pAb

Catalog No.: A13657

Basic Information

Observed MW

Calculated MW

51kDa

Category

Primary antibody

Applications

ELISA,IF/ICC

Cross-Reactivity

Human

Background

Enables tRNA methyltransferase activity. Predicted to be involved in tRNA wobble uridine modification. Predicted to be active in cytoplasm and nucleus.

Recommended Dilutions

IF/ICC 1:50 - 1:100

Immunogen Information

Gene ID

57604

Swiss Prot

Q9P272

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 110-330 of human KIAA1456 (NP_065895.2).

Synonyms

TRM9L; hTRM9L; C8orf79; KIAA1456

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.