

# IAH1 Rabbit pAb

Catalog No.: A13710

## Basic Information

### Observed MW

27kDa

### Calculated MW

28kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human, Mouse, Rat

## Background

Enables identical protein binding activity. Predicted to be involved in lipid catabolic process.

## Recommended Dilutions

WB 1:500 - 1:1000

## Immunogen Information

### Gene ID

285148

### Swiss Prot

Q2TAA2

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-248 of human IAH1 (NP\_001034702.1).

### Synonyms

IAH1

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

# Validation Data



Western blot analysis of various lysates, using IAHI1 Rabbit pAb (A13710) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 3s.