

A13774

Leader in Biomolecular Solutions for Life Science



# JAKMIP1 Rabbit pAb

Catalog No.: A13774

## Basic Information

### Observed MW

73kDa

### Calculated MW

73kDa

### Category

Mouse Monoclonal Antibody

### Applications

WB, ELISA

### Cross-Reactivity

Human, Mouse, Rat

## Background

Enables GABA receptor binding activity and RNA binding activity. Involved in cognition. Is extrinsic component of membrane. Part of ribonucleoprotein complex.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

152789

### Swiss Prot

Q96N16

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 337-626 of human JAKMIP1 (NP\_653321.1).

### Synonyms

JAMIP1; MARLIN1; Gababrbp; JAKMIP1

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using JAKMIP1 Rabbit pAb (A13774) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 90s.