

RBSN Rabbit pAb

Catalog No.: A13792

Basic Information

Observed MW

Refer to figures

Calculated MW

89kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Mouse

Background

This gene encodes a protein that belongs to the FYVE zinc finger family of proteins. The encoded protein interacts with Ras-related proteins that regulate membrane trafficking. A missense mutation in this gene is associated with a defect in the early endocytic pathway. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

64145

Swiss Prot

Q9H1K0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human RBSN (NP_071735.2).

Synonyms

ZFYVE20; Rabenosyn-5; RBSN

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.