

DYNLL2 Rabbit pAb

Catalog No.: A13889

Basic Information

Observed MW

13kDa

Calculated MW

10kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Mouse, Rat

Background

Predicted to enable dynein intermediate chain binding activity and dynein light intermediate chain binding activity. Predicted to be involved in cilium assembly. Located in 9+0 non-motile cilium and centrosome. Is active in glutamatergic synapse and postsynapse.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

140735

Swiss Prot

Q96FJ2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-89 of human DYNLL2 (NP_542408.1).

Synonyms

Dlc2; DNCL1B; RSPH22; DYNLL2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using DYNLL2 Rabbit pAb (A13889) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 90s.