

A13909

Leader in Biomolecular Solutions for Life Science



# SNAP23 Rabbit pAb

Catalog No.: A13909

## Basic Information

### Observed MW

23kDa

### Calculated MW

23kDa

### Category

Polyclonal Antibody

### Applications

WB,IHC-P,ELISA

### Cross-Reactivity

Human,Mouse,Rat

## Background

Specificity of vesicular transport is regulated, in part, by the interaction of a vesicle-associated membrane protein termed synaptobrevin/VAMP with a target compartment membrane protein termed syntaxin. These proteins, together with SNAP25 (synaptosome-associated protein of 25 kDa), form a complex which serves as a binding site for the general membrane fusion machinery. Synaptobrevin/VAMP and syntaxin are believed to be involved in vesicular transport in most, if not all cells, while SNAP25 is present almost exclusively in the brain, suggesting that a ubiquitously expressed homolog of SNAP25 exists to facilitate transport vesicle/target membrane fusion in other tissues. The protein encoded by this gene is structurally and functionally similar to SNAP25 and binds tightly to multiple syntaxins and synaptobrevins/VAMPs. It is an essential component of the high affinity receptor for the general membrane fusion machinery and is an important regulator of transport vesicle docking and fusion. Two alternative transcript variants encoding different protein isoforms have been described for this gene.

## Recommended Dilutions

WB	1:500 - 1:2000
IHC-P	1:100 - 1:200

## Immunogen Information

Gene ID	Swiss Prot
8773	000161

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human SNAP23 (NP\_003816.2).

### Synonyms

SNAP-23; SNAP23A; SNAP23B; HsT17016; SNAP23

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information


Source	Isotype	Purification
Rabbit	IgG	Affinity purification

### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.


## Validation Data




Western blot analysis of various lysates using SNAP23 Rabbit pAb (A13909) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 2min.



Immunohistochemistry analysis of SNAP23 in paraffin-embedded rat lymph node using SNAP23 Rabbit pAb (A13909) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of SNAP23 in paraffin-embedded human appendix using SNAP23 Rabbit pAb (A13909) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of SNAP23 in paraffin-embedded mouse spleen using SNAP23 Rabbit pAb (A13909) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.