

# ZNF544 Rabbit pAb

Catalog No.: A14120

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

82kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse, Rat

## Background

---

Predicted to enable DNA-binding transcription activator activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be active in nucleus.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

27300

### Swiss Prot

Q6NX49

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ZNF544 (NP\_055295.2).

### Synonyms

ZNF544

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.