

# ROBO4 Rabbit pAb

Catalog No.: A14134

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

107kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human, Mouse

## Background

---

Predicted to enable cell-cell adhesion mediator activity. Involved in angiogenesis and establishment of endothelial barrier. Located in extracellular exosome. Implicated in aortic valve disease 3.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

54538

### Swiss Prot

Q8WZ75

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 848-1007 of human ROBO4 (NP\_061928.4).

### Synonyms

MRB; AOVD3; ECSM4; ROBO4

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.