

# PHC3 Rabbit pAb

Catalog No.: A14151

## Basic Information

---

### Observed MW

### Calculated MW

106kDa

### Category

Primary antibody

### Applications

ELISA,IF/ICC

### Cross-Reactivity

Human

## Background

---

Predicted to enable chromatin binding activity and histone binding activity. Predicted to be involved in negative regulation of transcription, DNA-templated. Located in nucleoplasm. Part of PRC1 complex.

## Recommended Dilutions

---

IF/ICC 1:50 - 1:100

## Immunogen Information

---

### Gene ID

80012

### Swiss Prot

Q8NDX5

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human PHC3 (NP\_079223.3).

### Synonyms

EDR3; HPH3; PHC3

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Immunofluorescence analysis of U2OS cells using PHC3 antibody (A14151). Blue: DAPI for nuclear staining.