

A14164

Leader in Biomolecular Solutions for Life Science



## PLCXD2 Rabbit pAb

Catalog No.: A14164

### Basic Information

#### Observed MW

#### Calculated MW

35kDa

#### Category

Mouse Monoclonal Antibody

#### Applications

WB,IF/ICC,ELISA

#### Cross-Reactivity

Human,Mouse,Rat

### Background

Predicted to enable phosphoric diester hydrolase activity. Predicted to be involved in lipid catabolic process and signal transduction. Predicted to be located in nucleus.

### Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:100

### Immunogen Information

#### Gene ID

257068

#### Swiss Prot

Q0VAA5

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human PLCXD2 (NP\_001172035.1).

#### Synonyms

PLCXD2

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification


#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Immunofluorescence analysis of A549 cells  
using PLCXD2 antibody (A14164).