

A14290

Leader in Biomolecular Solutions for Life Science

ZSCAN31 Rabbit pAb

Catalog No.: A14290



Basic Information

Observed MW

47kDa

Calculated MW

47kDa

Category

Mouse Monoclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Rat

Background

This gene encodes a protein containing multiple C2H2-type zinc finger motifs. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

64288

Swiss Prot

Q96LW9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 202-406 of human ZSCAN31 (NP_112161.3).

Synonyms

ZNF323; ZNF310P; ZNF20-Lp; ZSCAN31

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of lysates from rat liver, using ZSCAN31 Rabbit pAb (A14290) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.