

RFTN2 Rabbit pAb

Catalog No.: A14340

Basic Information

Observed MW

Refer to Figures

Calculated MW

56kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to act upstream of or within dsRNA transport and response to exogenous dsRNA. Predicted to be located in plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

130132

Swiss Prot

Q52LD8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 412-501 of human RFTN2 (NP_653230.2).

Synonyms

C2orf11; Raftlin-2; RFTN2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.