

# RPUSD2 Rabbit pAb

Catalog No.: A14397

## Basic Information

---

### Observed MW

Refer to Figures

### Calculated MW

61kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse, Rat

## Background

---

Enables pseudouridine synthase activity. Involved in mRNA pseudouridine synthesis.

## Recommended Dilutions

---

**WB** 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

27079

### Swiss Prot

Q8IZ73

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 286-545 of human RPUSD2 (NP\_689473.1).

### Synonyms

PUS9; C18B11; C15orf19; RPUSD2

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.