

# TEKT2 Rabbit pAb

Catalog No.: A14398

## Basic Information

---

### Observed MW

Refer to Figures

### Calculated MW

50kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human, Mouse, Rat

## Background

---

This gene product belongs to the tektin family of proteins. Tektins comprise a family of filament-forming proteins that are coassembled with tubulins to form ciliary and flagellar microtubules. This gene is expressed in the testis and its protein is localized to the flagella of the sperms, indicating that it may play a role in spermatogenesis.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

27285

### Swiss Prot

Q9UIF3

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-430 of human TEKT2 (NP\_055281.2).

### Synonyms

TEKTB1; TEKIN-T; h-tektin-t; TEKT2

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.