

# PCNP Rabbit pAb

Catalog No.: A14413

## Basic Information

---

### Observed MW

19kDa

### Calculated MW

19kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human

## Background

---

Involved in proteasome-mediated ubiquitin-dependent protein catabolic process and protein ubiquitination. Located in nuclear body.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

57092

### Swiss Prot

Q8WW12

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human PCNP (NP\_065090.1).

### Synonyms

PCNP

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of various lysates using PCNP Rabbit pAb (A14413) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 90s.