

NLGN2 Rabbit pAb

Catalog No.: A14418

Basic Information

Observed MW

Refer to Figures

Calculated MW

91kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Mouse, Rat

Background

This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

57555

Swiss Prot

Q8NFZ4

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 15-180 of human NLGN2 (NP_065846.1).

Synonyms

NLGN2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.