

ZNF585A Rabbit pAb

Catalog No.: A14457

Basic Information

Observed MW

Refer to Figures

Calculated MW

88kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to enable DNA-binding transcription activator activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be active in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

199704

Swiss Prot

Q6P3V2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human ZNF585A (NP_001275729.1).

Synonyms

Zfp27; ZNF585A

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.