

ASF1A Rabbit pAb

Catalog No.: A14508

Basic Information

Observed MW

Refer to Figures

Calculated MW

23kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse

Background

This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

25842

Swiss Prot

Q9Y294

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 100-204 of human ASF1A (NP_054753.1).

Synonyms

CIA; CGI-98; HSPC146; ASF1A

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.