

A14517

Leader in Biomolecular Solutions for Life Science



LNX1 Rabbit pAb

Catalog No.: A14517

Basic Information

Observed MW

Refer to Figures

Calculated MW

81kDa

Category

Mouse Monoclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

This gene encodes a membrane-bound protein that is involved in signal transduction and protein interactions. The encoded product is an E3 ubiquitin-protein ligase, which mediates ubiquitination and subsequent proteasomal degradation of proteins containing phosphotyrosine binding (PTB) domains. This protein may play an important role in tumorigenesis. Alternatively spliced transcript variants encoding distinct isoforms have been described. A pseudogene, which is located on chromosome 17, has been identified for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

84708

Swiss Prot

Q8TBB1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 480-620 of human LNX1 (NP_001119800.1).

Synonyms

LNX; MPDZ; PDZRN2; LNX1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.