

# PWP1 Rabbit pAb

Catalog No.: A14579

## Basic Information

---

### Observed MW

Refer to Figures

### Calculated MW

56kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse

## Background

---

The protein encoded by this gene contains several WD-40 repeats and is found mostly in the nucleus. The expression and localization of this protein are cell cycle dependent. Expression of this gene is upregulated in pancreatic adenocarcinoma. Three transcript variants encoding different isoforms have been found for this gene.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

11137

### Swiss Prot

Q13610

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human PWP1 (NP\_008993.1).

### Synonyms

NCLB; IEF-SSP-9502; PWP1

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.