

A14587

Leader in Biomolecular Solutions for Life Science



DCUN1D1 Rabbit pAb

Catalog No.: A14587

Basic Information

Observed MW

30kDa

Calculated MW

30kDa

Category

Mouse Monoclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse

Background

Enables cullin family protein binding activity. Involved in positive regulation of protein neddylation and regulation of protein ubiquitination. Located in cytosol and nucleoplasm. Part of ubiquitin ligase complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

54165

Swiss Prot

Q96GG9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 15-70 of human DCUN1D1 (NP_065691.2).

Synonyms

RP42; SCRO; Tes3; DCNL1; SCCRO; DCUN1L1; DCUN1D1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using DCUN1D1 Rabbit pAb (A14587) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.