

A15018

Leader in Biomolecular Solutions for Life Science



ACAA1 Rabbit pAb

Catalog No.: A15018

Basic Information

Observed MW

Refer to figures

Calculated MW

44kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

This gene encodes an enzyme operative in the beta-oxidation system of the peroxisomes. Deficiency of this enzyme leads to pseudo-Zellweger syndrome. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

30

Swiss Prot

P09110

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 27-300 of human ACAA1 (NP_001598.1).

Synonyms

ACAA; THIO; PTHIO; Lnc-Myd88; ACAA1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.