

A15021

Leader in Biomolecular Solutions for Life Science



ACYP1 Rabbit pAb

Catalog No.: A15021

Basic Information

Observed MW

Calculated MW
11kDa

Category
Polyclonal Antibody

Applications
IF/ICC,ELISA

Cross-Reactivity
Human

Background

This gene is a member of the acylphosphatase family. The encoded protein is a small cytosolic enzyme that catalyzes the hydrolysis of the carboxyl-phosphate bond of acylphosphates. Two isoenzymes have been isolated and described based on their tissue localization: erythrocyte (common) type acylphosphatase encoded by this gene, and muscle type acylphosphatase. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID 97	Swiss Prot P07311
----------------------	-----------------------------

Immunogen
Recombinant fusion protein containing a sequence corresponding to amino acids 1-99 of human ACYP1 (NP_001098.1).

Synonyms
ACYPE; ACYP1

Contact


 | www.abclonal.com

Product Information

Source Rabbit	Isotype IgG	Purification Affinity purification
-------------------------	-----------------------	--

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of HeLa cells using ACYP1 Rabbit pAb (A15021) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (A5007) at 1:500 dilution. Blue: DAPI for nuclear staining.