

A15182

Leader in Biomolecular Solutions for Life Science



AEN Rabbit pAb

Catalog No.: A15182

Basic Information

Observed MW

Refer to figures

Calculated MW

36kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Enables exonuclease activity. Involved in intrinsic apoptotic signaling pathway in response to DNA damage by p53 class mediator and response to ionizing radiation. Located in nuclear membrane; nucleolus; and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

64782

Swiss Prot

Q8WTP8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human AEN (NP_073604.3).

Synonyms

ISG20L1; pp12744; AEN

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.