

A15901

Leader in Biomolecular Solutions for Life Science



ITFG1 Rabbit pAb

Catalog No.: A15901

Basic Information

Observed MW

80kDa

Calculated MW

68kDa

Category

Mouse Monoclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

81533

Swiss Prot

Q8TB96

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 313-612 of human ITFG1 (NP_110417.2).

Synonyms

TIP; CDA08; LNKN-1; 2310047C21Rik; ITFG1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using ITFG1 Rabbit pAb (A15901) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 3s.