

A15903

Leader in Biomolecular Solutions for Life Science



Defb26 Rabbit pAb

Catalog No.: A15903

Basic Information

Observed MW

Calculated MW

14kDa

Category

Mouse Monoclonal Antibody

Applications

IF/ICC,ELISA

Cross-Reactivity

Human

Background

Predicted to be involved in defense response to bacterium and innate immune response. Predicted to be located in extracellular region. Orthologous to human DEFB125 (defensin beta 125).

Recommended Dilutions

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID

641629

Swiss Prot

Q32ZG6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 22-128 of rat Defb26 (NP_001032595.1).

Synonyms

DEFB126; C20orf8; DEFB-26; DEFB26; HBD26; bA530N10.1; hBD-26; Defb26

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of U2OS cells using Defb26 antibody (A15903) at dilution of 1:100. Blue: DAPI for nuclear staining.