

A15903

Leader in Biomolecular Solutions for Life Science



## Defb26 Rabbit pAb

Catalog No.: A15903

### Basic Information

#### Observed MW

**Calculated MW**  
14kDa

**Category**  
Polyclonal Antibody

**Applications**  
IF/ICC,ELISA

**Cross-Reactivity**  
Human

### Background

Predicted to be involved in defense response to bacterium and innate immune response. Predicted to be located in extracellular region. Orthologous to human DEFB125 (defensin beta 125).

### Recommended Dilutions

IF/ICC 1:50 - 1:200

### Immunogen Information

Gene ID	Swiss Prot
641629	Q32ZG6

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 22-128 of rat Defb26 (NP\_001032595.1).

#### Synonyms

DEFB126; C20orf8; DEFB-26; DEFB26; HBD26; bA530N10.1; hBD-26; Defb26

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information


Source	Isotype	Purification
Rabbit	IgG	Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Immunofluorescence analysis of U2OS cells using Defb26 antibody (A15903) at dilution of 1:100. Blue: DAPI for nuclear staining.