

A15913

Leader in Biomolecular Solutions for Life Science



ACBD6 Rabbit pAb

Catalog No.: A15913

Basic Information

Observed MW

Refer to figures

Calculated MW

31kDa

Category

Mouse Monoclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to enable fatty-acyl-CoA binding activity and lipid binding activity. Predicted to be located in cytosol.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

84320

Swiss Prot

Q9BR61

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-282 of human ACBD6 (NP_115736.1).

Synonyms

ACBD6

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.