

SRRM2 Rabbit pAb

Catalog No.: A15993

Basic Information

Observed MW

Refer to figures

Calculated MW

300kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

Enables C2H2 zinc finger domain binding activity and protein N-terminus binding activity. Involved in mRNA splicing, via spliceosome. Located in Cajal body and nuclear speck. Part of U2-type catalytic step 2 spliceosome and U2-type pre-catalytic spliceosome. Biomarker of Parkinson's disease.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

23524

Swiss Prot

Q9UQ35

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SRRM2 (NP_057417.3).

Synonyms

CWF21; Cwc21; 300-KD; SRL300; SRm300; HSPC075; SRRM2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.