

A16025

Leader in Biomolecular Solutions for Life Science



# FRMD3 Rabbit pAb

Catalog No.: A16025

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

69kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human

## Background

---

The protein encoded by this gene is a single pass membrane protein primarily found in ovaries. A similar protein in erythrocytes helps determine the shape of red blood cells, but the function of the encoded protein has not been determined. There is some evidence that this is a tumor suppressor gene, and there is also evidence linking defects in this gene to susceptibility to diabetic nephropathy in type 1 diabetes. Several transcript variants encoding different isoforms have been found for this gene.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

257019

### Swiss Prot

A2A2Y4

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of human FRMD3 (NP\_777598.3).

### Synonyms

4.10; P410; EPB41LO; EPB41L40; FRMD3

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.