

A17385

Leader in Biomolecular Solutions for Life Science



CSNK1A1L Rabbit pAb

Catalog No.: A17385

Basic Information

Observed MW

Refer to figures

Calculated MW

39kDa

Category

Polyclonal Antibody

Applications

WB,IF/ICC,IP

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to enable protein serine/threonine kinase activity. Predicted to be involved in negative regulation of canonical Wnt signaling pathway and peptidyl-serine phosphorylation. Located in cytosol.

Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200
IP	1:50 - 1:200

Immunogen Information

Gene ID	Swiss Prot
122011	Q8N752

Immunogen

Recombinant Protein of human CSNK1A1L.

Synonyms

CK1; CSNK1A1L

Contact

 www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.