

A17394

Leader in Biomolecular Solutions for Life Science



PPP1R13B Rabbit pAb

Catalog No.: A17394

Basic Information

Observed MW

Refer to figures

Calculated MW

120kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a member of the ASPP (apoptosis-stimulating protein of p53) family of p53 interacting proteins. The protein contains four ankyrin repeats and an SH3 domain involved in protein-protein interactions. ASPP proteins are required for the induction of apoptosis by p53-family proteins. They promote DNA binding and transactivation of p53-family proteins on the promoters of proapoptotic genes. Expression of this gene is regulated by the E2F transcription factor.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

23368

Swiss Prot

Q96KQ4

Immunogen

A synthetic peptide of human PPP1R13B.

Synonyms

p85; ASPP1; p53BP2-like; PPP1R13B

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.