

A17751

Leader in Biomolecular Solutions for Life Science



ABHD4 Rabbit pAb

Catalog No.: A17751

Basic Information

Observed MW

39kDa

Calculated MW

39kDa

Category

Mouse Monoclonal Antibody

Applications

WB,IF/ICC,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to enable lysophosphatidic acid acyltransferase activity and lysophospholipase activity. Predicted to be involved in N-acylphosphatidylethanolamine metabolic process; lipid homeostasis; and phosphatidic acid biosynthetic process. Predicted to act upstream of or within N-acylethanolamine metabolic process. Predicted to be located in endoplasmic reticulum membrane. Predicted to be active in lipid droplet and mitochondrion.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID

63874

Swiss Prot

Q8TB40

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 223-342 of human ABHD4 (NP_071343.2).

Synonyms

ABH4; ABHD4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of lysates from Rat testis, using ABHD4 Rabbit pAb (A17751) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.



Immunofluorescence analysis of L929 cells using ABHD4 Rabbit pAb (A17751) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.