

ADAD2 Rabbit pAb

Catalog No.: A17840

Basic Information

Observed MW

Refer to figures

Calculated MW

62kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to enable double-stranded RNA adenosine deaminase activity; double-stranded RNA binding activity; and tRNA-specific adenosine deaminase activity. Predicted to be involved in RNA processing and adenosine to inosine editing. Predicted to be active in cytoplasm and nucleolus.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

161931

Swiss Prot

Q8NCV1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ADAD2 (NP_001138872.1).

Synonyms

TENRL; ADAD2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.