

A17882

Leader in Biomolecular Solutions for Life Science

## ABCB5 Rabbit pAb

Catalog No.: A17882



### Basic Information

---

**Observed MW**

Refer to figures

**Calculated MW**

139kDa

**Category**

Polyclonal Antibody

**Applications**

ELISA

**Cross-Reactivity**

Human, Mouse

### Background

---

ABCB5 belongs to the ATP-binding cassette (ABC) transporter superfamily of integral membrane proteins. These proteins participate in ATP-dependent transmembrane transport of structurally diverse molecules ranging from small ions, sugars, and peptides to more complex organic molecules (Chen et al., 2005 [PubMed 15760339]).

### Recommended Dilutions

---

### Immunogen Information

---

**Gene ID**

340273

**Swiss Prot**

Q2M3G0

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 120-230 of human ABCB5 (NP\_848654.3).

**Synonyms**

ABCB5beta; EST422562; ABCB5alpha; ABCB5

### Contact

---



[www.abclonal.com](http://www.abclonal.com)

### Product Information

---

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.