

A17986

Leader in Biomolecular Solutions for Life Science



MonoMethyl-eEF1A-K55 Rabbit pAb

Catalog No.: A17986

Basic Information

Observed MW

Calculated MW

Category

Mouse Monoclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID **Swiss Prot**

Immunogen

A synthetic monomethylated peptide around K55 of human eEF1A (NP_001393.1).

Synonyms

Contact



www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit		

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.